

VISUAL URINALYSIS STRIPS 11 PARAMETERS - professional strips



Product Code	24077
Unit of sale	box of 100 pcs
Minimum order	1
Type	IVD - Professional Use
Class	

RDM (NSIS)	2085574
CND	W0101060204

Description

11 PARAMETERS STRIPS - tube of 100 strips
PROFESSIONAL USE
Professionals strips for testing and analysis of urine specimen with visual results. Strips can also be examined with Gima urine analyzer (24046).
- anti-VC interference ability;
- global check/screening for general health condition;
- rapid results: from 30 to 120 seconds;
- high quality colours chart;
- cut-off value measures for each parameter.
Multilingual instructions: GB, FR, IT, ES, PT, SE, HU, GR.

Technical Specifications

Parameters

Glucose: Testing and reading time 30 s, Cut-off value 2.8 ~ 5.5 mmol/L
Ascorbic Acid: Testing and reading time 60 s, Cut-off value 0.3 ~ 0.6 mmol/L
Ketone: Testing and reading time 40 s, Cut-off value 0.5 ~ 1.0 mmol/L
Protein: Testing and reading time 60 s, Cut-off value 0.15 ~ 0.3 g/L
pH-Value: Testing and reading time 60 s, Cut-off value 0.5
Blood: Testing and reading time 60 s, Cut-off value 5 - 15 Ery/ μ L
Nitrite: Testing and reading time 60 s, Cut-off value 13 ~ 22 μ mol/L
Leucocytes: Testing and reading time 120 s, Cut-off value 5 - 15 Leuko/ μ L
Specific Gravity: Testing and reading time 45 s, Cut-off value 0.005
Bilirubin: Testing and reading time 30 s, Cut-off value 8.6 - 17 μ mol/L
Urobilinogen: Testing and reading time 60 s, Cut-off value 3.3 - 16 μ mol/L