

CONTOUR HEAD IMMOBILIZER - black/yellow



Product Code	34079
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2327938
CND	Y060312
EAN/UPC	8023279340792

Description

HEAD IMMOBILIZERS

CONTOUR - adult
Colour: yellow
Weight (kg): 1.5
Material: moulded polyurethane and PVC

Adult and child models of head immobilizers, different in shape, material and colour. All are composed of two head blocks and two straps to prevent any head movement of the injured person during transport. Both sides have an ear hole. They can be used with scooping and PE stretchers.

Radiolucent material, suitable for x-ray chamber

Multilingual manual: GB, FR, IT, ES, PT, DE, NL, DK, GR, BG, RO