

CARRYING COT 190x75 cm



Product Code	34090
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2534389
CND	V08050199
EAN/UPC	8023279340907
UMDNS	13818
PARAF	903163424
GTIN	08023279340907
GMDN	65516

Description

Made in polyester, PVC coated.
Anti-putridity, bacterial and fungus proof, waterproof, easy to wash and disinfect, X-ray translucent material.

Size 190 x 75 cm.
Loading capacity: 150 kg.
Certified EN1865-2010 + A1 2015
Supplied in a carrying bag.

Multilanguage manual: GB, FR, IT, ES, PT, DE.

With 8 handles for a secure patient transport.