



Product Code	41001
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	1234580
CND	T020401
EAN/UPC	8023279410013

Description

100% TEXTURIZED POLYESTER DRAPES

- POLYESTER DRAPE 150x150 cm - light blue

Studied for replacing the classic cotton drape in operating rooms.

Advantages:

- it does not have particles release of dust or cotton fibers that could cause infections during surgery

Various usability on the field:

- as a sheet under the patient
- as a towel
- to tamponage, considering its high absorbency
- as a pre or post surgery cover, thanks to the capability of keeping a good body temperature