



Product Code	34107
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2210087
CND	V0899
EAN/UPC	8023279341072

Description

GIMA STRETCHER 4 - foldable in 4

Made of strong alloy and watertight material. Light weight, easy to carry, foldable in 4. Length foldable in 2 parts, width foldable in 2 parts.

Supplied with carrying bag.

Technical Specifications

- Open size: 221 x 53 x 16 cm
- Folded size: 104 x 19 x 12 cm
- Weight: 5 kg
- Load: 150 kg