



Product Code	24111
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No
RDM (NSIS)	1832273
CND	W0101050205
EAN/UPC	8023279241112
PARAF	978866010

Description

GIMA GLUCOSE MONITOR mg/dL

Glucose monitor kit

Box and manual languages: IT, DE, GR, Arabic

Accurate and reliable glucose monitors for self testing and professional use.

Extensive clinical studies demonstrate the accuracy of Gima glucose monitoring system, powered by Acon, with fresh capillary, venous, neonatal blood samples.

Comply with ISO15197:2013/EN ISO15197:2015.

Kit includes 10 strips, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual.

Technical Specifications

Technology: Biosensor/Electrochemical, Glucose dehydrogenase (GDH-FAD)

Calibration method: No coding

Result calibration: Plasma-equivalent

Test time: 5 seconds

Sample size: 0.6µL

Sample type: Fresh whole blood (capillary, venous, neonatal)

Test sites: Fingertip, forearm, palm

Hematocrit range: 10-70%

Memory storage: 1,000 results with date and time

Test averaging: 7, 14, 30, 60, 90-day averages

Data transfer: USB cable - optional

Control solution: 3 levels

Automatic shutoff: 2 minutes after last action

Battery: 2x CR 2032 3.0V coin cell batteries

Battery life: 3,000 measurements

Other functions: Automatic detection of insufficient sample / 2nd sample application allowed / Sample perturbation detection