

# STAY UPS - XL - medium compression - beige



|               |                |
|---------------|----------------|
| Product Code  | 41127          |
| Unit of sale  | 1 pair         |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2636026       |
| CND        | M030405       |
| EAN/UPC    | 8033201400463 |
| PARAF      | 925564229     |

## Description

Stay ups with gradual compression, silicone lace top, honeycomb mesh

- Medium compression Stay Ups
- Colour: Beige
- Size: X-Large

Stay Ups designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics.  
This line offers two different compression levels (medium and strong) and is available in beige or black.

- Medium compression level: 70 den - 10-14 mm/Hg at the ankle Ideal for those who would like to initiate a prevention programme, but are mainly interested in the comfort and well-being offered by elastic compression stockings.

Composition: 84% polyamide, 16% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.