



Product Code	41137
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2636027
CND	M030405
EAN/UPC	8033201400777
PARAF	925565285

## Description

Stay ups with gradual compression, silicone lace top, honeycomb mesh

- Strong compression Stay Ups
- Colour: Beige
- Size: Medium

Stay Ups designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics.

This line offers two different compression levels (medium and strong) and is available in beige or black.

- Strong compression level: 140 den - 15-21 mm/Hg at the ankle

Ideal for those suffering from the initial signs of ankle swelling and discomfort. This level of compression is recommended for initial onset of varicose veins. The product's compression level supports and massages the leg, providing effective relief throughout the day. While this product line boasts an elegant design, particular attention has also been paid to maximum elasticity, which is ensured through the use of a honeycomb fabric (net fabric).

Doctor consult is recommended before use.

Composition: 79% polyamide, 21% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.