

# GIMA GLUCOSE MONITOR KIT mg/dL with Bluetooth - GB,FR,IT,ES,PT,DE



Product Code	24114
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

UK-REP	Yes
CH-REP	No

RDM (NSIS)	1853862
CND	W0101050205
EAN/UPC	8023279241143
PARAF	980907051

## Description

GIMA GLUCOSE MONITOR mg/dL

**Glucose monitor kit with Bluetooth**  
**Box and manual languages: GB, FR, IT, ES, PT, DE**

Accurate and reliable glucose monitors for self testing and professional use. Extensive clinical studies demonstrate the accuracy of Gima glucose monitoring system, powered by Acon, with fresh capillary, venous, neonatal blood samples. Comply with ISO15197:2013/EN ISO15197:2015.

Kit includes 10 strips, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual.

**Bluetooth version**  
Fully integrated Bluetooth allows for data transmission to data management platforms.

## Technical Specifications

Technology: Biosensor/Electrochemical, Glucose dehydrogenase (GDH-FAD)  
Calibration method: No coding  
Result calibration: Plasma-equivalent  
Test time: 5 seconds  
Sample size: 0.6µL  
Sample type: Fresh whole blood (capillary, venous, neonatal)  
Test sites: Fingertip, forearm, palm  
Hematocrit range: 10-70%  
Memory storage: 1,000 results with date and time  
Test averaging: 7, 14, 30, 60, 90-day averages  
Data transfer: USB cable - optional  
Control solution: 3 levels  
Automatic shutoff: 2 minutes after last action  
Battery: 2x CR 2032 3.0V coin cell batteries  
Battery life: 3,000 measurements  
Other functions: Automatic detection of insufficient sample / 2nd sample application allowed / Sample perturbation detection