

# BODY BAG vinyl covered nylon - load 150kg



Product Code	34117
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279341171
---------	---------------

PARAF	977813967
-------	-----------

## Description

BODY BAG - vinyl - blue

Made of vinyl covered nylon conform to standard EN 1865. Provided with a long slide fastener and 6 handles.

Size: 210x90 cm.

Max load: 150 kg.