

PANTYHOSES - XXL - strong compression - beige



Product Code	41174
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2636030
CND	M030405
EAN/UPC	8033201400883
PARAF	925565400

Description

Pantyhose with gradual compression, reinforced body, honeycomb mesh

- Strong compression Pantyhoses
- Colour: Beige
- Size: XX-Large

Pantyhoses designed for dynamic women who need the preventive action of elastic compression but don't want to compromise on style.

This elegant product line features a true heel, an anti-shrink cuff and extremely soft fabrics.
This line offers two different compression levels (medium and strong) and is available in beige or black.

- Strong compression level: 140 den - 15-21 mm/Hg at the ankle
Ideal for those suffering from the initial signs of ankle swelling and discomfort. This level of compression is recommended for initial onset of varicose veins. The product's compression level supports and massages the leg, providing effective relief throughout the day. While this product line boasts an elegant design, particular attention has also been paid to maximum elasticity, which is ensured through the use of a honeycomb fabric (net fabric).

Doctor consult is recommended before use.

Composition: 79% polyamide, 21% elastane.

Multilanguage box and instructions: GB, FR, IT, ES, PT, DE, GR, HU

Made in Italy.