

FULL LENGTH STOCKINGS length 70-80 - X-large



Product Code	41333
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2636044
CND	M030405
EAN/UPC	8033201404973
PARAF	986897698

Description

Antiembolism full length stocking with inspection hole adjustable top, ambidextrous

- Size: X-Large
- Length: 70-80 cm
- Compression mm/Hg 18-20 (at the ankle)

A line of stockings dedicated to prophylaxis of thromboembolism before, during and after surgery.

Thanks to their descending graduated compression, these stockings represent an excellent means of preventing increased blood viscosity, which can lead to the formation of thrombus.

In order to guarantee a good fit and properly graduated compression, a great deal of attention has been dedicated to ensuring the high elasticity of the fabrics.

- Indications:
 - a valid aid in preventing thromboembolisms linked to surgical interventions.
 - treatments associated with post-traumatic therapies.
 - bruises, sprains and hydrostatic edemas

- Contraindications:
 - localized dermatitis, erythema, necrosis
 - stripping, pulmonary edema, grafts
 - vascular disorders of cardiac origin
 - legs of larger circumference

- Features:
 - decreasing compression from the ankle to the thigh
 - true anatomic heel
 - extremely soft and elastic fabric, hypoallergenic (latex free)
 - inspection opening on top of the foot

Composition: 82% polyamide, 18% elastane.

Multilanguage bag and instructions: GB, FR, IT, ES, PT, DE, GR, HU.

Made in Italy.