



Product Code	24147
Unit of sale	1 pc
Minimum order	1
Type	IVD - Self Test
Class	

CND	W0299
EAN/UPC	8032937231259

Description

MULTICAREIN - Italian

Supplied with nylon bag 5 cholesterol strips, 5 safety lancets, rapid guide (GB, FR, IT, ES, PT, GR, PL, HU) and IT manual.

Multiparameter diagnostic device for the periodical self-monitoring of risk factors: cholesterol, triglycerides and glucose. Portable device for easy test procedure and quick results.

Accurate with Amperometric Technology for glucose and Reflectometric Technology for cholesterol and tryglicerides.

Switchable mg/dL - mmol/L

Multilingual box: GB, FR, IT, ES, GR, PT, PL, HU, RU and selectable internal software: GB, IT.

Technical Specifications

- Glucose electro-chemical
Measuring range: 10-600 mg/dL - 0.6-33.3 mmol/L
Measuring time: 5"
Blood sample: 0.5 µl
- Cholesterol
Measuring range: 130-400 mg/dL - 3-10.2 mmol/L
Measuring time: 30"-60"
Blood sample: 10 µl
- Triglycerides
Measuring range: 50-500 mg dL - 0.56 - 5.6 mmol/L
Measuring time: 30"-60"
Blood sample: 10 µl
- Memory: 500 test results with date and time
- Display: LCD
- Batteries: 2 lithium batteries Cr 2032 3 Volts (32212) suitable for 1 year or 1,000 measurements
- Size - weight: 97x49x20.5 mm - 65 g (battery included)