

# SINGLE ELECTRIC BREAST PUMP



|               |                |
|---------------|----------------|
| Product Code  | 41701          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | II A           |

|            |               |
|------------|---------------|
| UK-REP     | No            |
| CH-REP     | No            |
| RDM (NSIS) | 2001713       |
| CND        | Z12080303     |
| EAN/UPC    | 6970392218908 |
| UMDNS      | 10485         |
| PARAF      | 984594883     |

## Description

### SINGLE ELECTRIC BREAST PUMP

Compact and lightweight breast pumps for efficient and discreet pumping. Easy to assemble, use and clean. Supplied with breast shield, nipple with cover, milk bottle with lid, multilingual box and manual (GB, FR, IT, ES, PT, DE).

All parts that touch breast milk are made without BPA, phthalate free.  
2 phase design, easy for stimulation and expression.  
10 vacuum levels  
Soft and comfortable breast shield  
Graduated bottle (ml/Oz)

## Technical Specifications

Vacuum Stimulate: - 60-130 mmHg  
Vacuum Expression: - 60-250 mmHg  
Weight: 265 g  
Size: 129 x 129 x 55 mm  
Operation: 4 AA batteries not included or AC adapter (optional)