

# "GIMA 1" EMERGENCY DURABLE CASE



Product Code	34176
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279341768
PARAF	977814058
GTIN	08023279341768

## Description

### EMERGENCY KITS GIMA 1

Durable case - orange (27222)  
Size: 55 x 43 x h 21 cm

Products included in kits:

- Emergency cover (34085)
- Universal scissors - 16.5 cm (20571)
- Facemask n. 2 - child/small (34222)
- Silicone resuscitator bag with mask n 4 - adult (34245)
- Teeth screw driver (34273)
- Tongue holding forceps (34274)
- Oxygen tube - 120 cm (34276)
- Mixed Guedel airways - box of 10 (34439)
- Oxygen therapy mask - adult (34166)
- 0.5 l oxygen cylinder (empty) with reducer UNI std (34501)
- Jotarap dual head stethoscope - black (32580)
- London sphygmomanometer (32725)
- Laryngoscope set 3 Mc-Intosh blades 2-3-4 (34302)
- Hand suction unit (28125)
- Delta pen light (25629)
- Fast tourniquet - green (25728)
- Magill forceps - 20 cm (34002)
- Gima digital thermometer °C (25560)
- Disinfectant wipes - tube 100 wipes (36625)
- Emocontrol elastic bandage (34782)
- Gauze bandage 3.5 m x 7 cm - box of 10 (34841)
- Sterile adhesive dressing 10 x 10 cm - box of 5 (34819)
- Instant ice - bag 3 pcs. (34111)
- 3 E3 padded pouches with window red (27188), blue (27189), green (27191)
- 2 E4 pouches with window red (27192), yellow (27194)