

# SPEED-1 FIRST AID CASE with cylinder UNI - autoclavable resuscitator



Product Code	34181
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279341812
PARAF	903163917
GTIN	08023279341812

## Description

First aid kits for reanimation, very useful for medical and dental rooms, sport clubs, factories, ambulances...  
Case made of anti-shock plastic.  
Size: 43 x 32 x 11 cm.  
Made in Italy.

## Standard accessories

- Silicone resuscitation bag - autoclavable
- Silicone Face mask n° 4
- Oxygen cylinder 0.5 l (empty). UNI standard
- High pressure manometer
- Guedel airways (small + medium + large)
- Helicoidal mouthwedge + Tongue holding forceps
- Atoxing tubing for connection to oxygen cylinder
- Plastic case