

SPEED-3 FIRST AID CASE without cylinder



Product Code	34183
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279341836
PARAF	903163968
GTIN	08023279341836

Description

First aid kits for reanimation, very useful for medical and dental rooms, sport clubs, factories, ambulances. Case made of anti-shock plastic.
Size: 43 x 32 x 11 cm.
Made in Italy.

Standard accessories

- Silicone resuscitation bag - autoclavable
- Silicone Face mask n° 4
- Guedel airways (small + medium + large)
- Helicoidal mouthwedge + Tongue holding forceps
- Atoxig tubing for connection to oxygen cylinder
- Plastic case