

SPLINTER REMOVAL KIT



Product Code	34190
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	Kit of medical devices
------	------------------------

Class	
-------	--

RDM (NSIS)	137689
------------	--------

EAN/UPC	8034028010347
---------	---------------

PARAF	903164034
-------	-----------

Description

Includes everything to remove splinters and for first aid medications.
All contents are in a polypropylene box size 235 x 190 x 43 mm.

Standard accessories

- 1 Professional removal pencil
- 1 Mirror
- 1 Magnifying glass
- 1 Eye wash solution 100 ml with an eye glass
- 1 Stainless steel tweezers 11.5 cm
- 2 sterile gauze pads 18x40 cm