

MULTIRED LARGE KIT - plastic case 3



Product Code	34193
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
RDM (NSIS)	136318
EAN/UPC	8056389237720

Description

MULTIRED LARGE KIT - plastic case 3

Size: 443x338x147 mm

Made of strong tripolymer plastic which is shiny smooth and easy to clean.
Neopren seal fitted to make it dustproof.
A wall bracket secures it on any house wall or booth.
Internal separators and safe closure with 2 swivel catches.

Weight: 3.3 (empty), 8.4 kg (filled).

Standard accessories

- 5 pairs of sterile gloves
- 1 safety mask with eye protection visor
- 3 sterile saline solutions - 500 ml
- 2 Iodopovidone disinfectant solutions 10% - 500 ml
- 10 sterile gauzes in sachets 10x10 cm
- 2 sterile gauzes in sachets 18x40 cm
- 2 sterile bandages 40x60 cm
- 1 pack of cotton
- 2 sterile tweezers
- 1 pack of net elastic bandage
- 2 adhesive plasters 5 m x 2.5 cm
- 2 Plastosan 20 assorted dressings
- 1 pair of lister scissors 14.5 cm
- 3 tourniquets
- 2 disposable instant ICE packs
- 2 sanitary bags
- 1 clinical digital thermometer
- 1 sphygmomanometer with stethoscope
- 2 gauze bandages 10 cm x 3.5 m
- 1 isothermic blanket gold/silver 160x210 cm
- 1 triangular TNT bandage 136x96x96 cm
- 1 pack of 8 sachet soap/disinfectant

2 elastic bandages 4 m x 6 cm
1 sterile compressive bandage 80x100
1 Pinocchio mouth opener
2 burn gel sachets 3.5 g