

LARGE KIT - metal cabinet



Product Code	34199
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

RDM (NSIS)	142028
EAN/UPC	8056389237904
PARAF	922540582

Description

- Large kit - metal cabinet
- Size: 460 x 300 x 140 mm
 - Weight: 7.5 kg

Standard accessories

- 2 Bottle antiseptic solution 500 ml
- 5 Pair protective gloves, sterile
- 2 Box of 10 assorted plasters
- 2 Adhesive plaster roll 2.5 cm x 5 m
- 1 Tubular net gauze
- 10 Pack of sterile compress 10 x 10 cm
- 2 Pack of sterile compress 18 x 40 cm
- 1 Pack of cotton wool
- 1 Mask with face shield
- 1 Pair of scissors 14.5 cm
- 2 Sterile forceps
- 3 Haemostatic lacet
- 1 Sphygmomanometer + Stethoscope
- 2 Sterile blanket cm 40 x 60
- 2 Pack of instant ice, single use
- 1 Digital thermometer
- 2 Trash bag
- 3 Irrigating solution 500 ml