

FABRIC ARM SLING - S - blue



Product Code	42000
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2737106
CND	Y090612
EAN/UPC	8023279420005
PARAF	950080061

Description

- Arm sling-fabric - S
- Length: 25-28 cm

Arm sling designed for arm rest and immobilization.
Ideal for immobilization after rotator cuff repair or shoulder girdle fracture, shoulder dislocation and subluxation.
Breathable and durable material with adjustable padded shoulder strap.
Envelope style with thumb support.
Machine washable.
Fits left and right.
Dark blue colour.