

AMBU MARK IV RESUSCITATOR - adult - with reservoir
 bag 1,500 ml



Product Code	34201
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	Yes
EAN/UPC	5707480019434

Description

AMBU® MARK IV RESUSCITATOR - adult - with reservoir

Designed in accordance with the original double wall principle, the new 100% latex-free Ambu bags are more robust and provide users with exceptional visual and tactile lung condition feedback during manual resuscitation, making operation safer and more secure. Built-in pressure limitation prevents compression if there is any resistance in the patient's airways.

Suitable for adults and children over 15 kg (3 years).

Supplied with face mask N°5, reservoir bag 1500 ml and single shutter valve.

Autoclavable at 134°C.
 Maximum volume 1300 ml.