



Product Code	42011
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2472395
CND	Y090612
EAN/UPC	8023279420111
PARAF	950080111

Description

UNIVERSAL ELBOW BRACE WITH SLING - left

Elbow brace designed for post-operative immobilization of the elbow in 15° increments and limitation of ROM (range of motion) by adjustable dial lock hinges.

Ideal for:

- conservative treatment of elbow dislocation or luxation
- stable or internally fixed fractures of the distal humerus or proximal radius or ulna.

Features:

- length adjustment (34-43 cm), one size fits all
- extension limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°
- flexion limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- immobilization limitation at: 0°, 15°, 30°, 45°, 60°, 75°, 90°
- breathable laminated soft foam material offers air permeability and moisture vapor transport

Supplied with shoulder sling for extra support and relief of arm.