



Product Code	42021
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2737096
CND	Y090618
EAN/UPC	8023279420210
PARAF	950080150

## Description

- Knee support - M
- Thigh circumf. cm: 47-53.5
- Calf circumf. cm: 37-43

Knee support indicated for:

- mild to moderate instabilities of knee joint and strain of the collateral ligaments
- meniscus injuries
- chronic inflammation after operation or trauma

It provides collateral stabilization of knee joint preventing valgus or varus deviation and hyperextension.

Closed cell neoprene ensures compressive and proprioceptive effect.

Blood circulation increases through warm effect.

Pull-on version with soft patellar ring that gives stabilization and prevents patellar displacement; soft-cushioned felt forms the contour of knee offering comfort; upper and lower coronal fasteners offer further fit.