



Product Code	42026
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2737096
CND	Y090618
EAN/UPC	8023279420265
PARAF	950080123

Description

UNIVERSAL POST OPERATION KNEE BRACE

Knee brace ideal for post operative immobilization or ROM control of ligament and fracture-related surgeries, injured or operated collateral or cruciate ligaments.

It provides excellent collateral stabilization via long side bars preventing valgus or varus deviation.

Adjustable ROM by setting hinges.

Features:

- "One push button" can be released for length adjustment
- extension limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- flexion limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
- immobilization limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- length adjustments from 48 to 70 cm
- transcendent overall adjustability
- simply press and rotate settings

Circumference: lower leg 48 cm, thigh 70 cm.

Can be worn over pants or on skin.