



Product Code	34204
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

GTIN	05707480035199
------	----------------

Description

Designed in accordance with the original double wall principle, the new 100% latex-free Ambu bags are more robust and provide users with exceptional visual and tactile lung condition feedback during manual resuscitation, making operation safer and more secure. Built-in pressure limitation prevents compression if there is any resistance in the patient's airways.

Autoclavable at 134°C.

Suitable for neonates, infants and children up to 20 kg (approx 4-5 years).

Supplied with face mask N°0A, oxygen reservoir tube 250 mm and pressure limiting valve 4.0 kPa (40 cm H₂O).
Maximum volume 300 ml.