



| | |
|---------------|----------------|
| Product Code | 42050 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|------------|---------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2472349 |
| CND | Y090618 |
| EAN/UPC | 8023279420500 |
| PARAF | 950080186 |

Description

Ankle brace with Q-Fit closure system for excellent fit, convenience and durability.

- Ankle brace - S

Indicated for:

- mild ankle distortion
- chronic ankle instability
- prophylactic and therapeutic use in sport activities

Made of durable and tough nylon fabric with bilateral metal stays for additional support. Extra nonelastic cross-stirrup straps provide maximal stability.