

WRIST BRACE - left - XL



Product Code	42063
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2472383
CND	Y090612
EAN/UPC	8023279420630
PARAF	950080251

Description

Q-fit wrist brace with reel locker designed for immobilization of the wrist.

- Wrist brace - XL
- Left

Ideal for:

- carpal tunnel syndrome
- TFCC (triangular fibrocartilage complex) injury
- radial nerve paralysis.

Q-fit dial closure system guarantees excellent fit, convenience and durability.
Comfortable, breathable skin-friendly material with anatomically-contoured aluminium palmar splint.