



Product Code	42067
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2472383
CND	Y090612
EAN/UPC	8023279420678
PARAF	950080287

## Description

Q-fit wrist brace with reel locker designed for immobilization of the wrist.

- Wrist brace - L
- Right

Ideal for:

- carpal tunnel syndrome
- TFCC (triangular fibrocartilage complex) injury
- radial nerve paralysis.

Q-fit dial closure system guarantees excellent fit, convenience and durability.

Comfortable, breathable skin-friendly material with anatomically-contoured aluminium palmar splint.