



Product Code	42102
Unit of sale	box of 20 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2472370
CND	Y090699
EAN/UPC	8023279421026

Description

FINGER SPLINT

- Size: 3

Available in 5 different sizes in clear colour.

Made of high quality rigid plastic.

With holes for air circulation and cooling.

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, PL, DK, SE, FI, CZ, SK, HU, HR, SI, RO, BG, GR, LT, LV, EE, Arabic.