



Product Code	34236
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No
GTIN	05707480021543

Description

AMBU SPUR II PEDIATRIC - with face mask n°1 and reservoir bag

Ambu SPUR II provides users with exceptional tactile and visual feedback during resuscitation.

The bag is highly responsive, with minimal mechanical resistance. Furthermore the characteristic design provides optimum stroke volume with perfect recoil.

Range includes infant, pediatric, and adult sizes.

All Ambu SPUR II resuscitators come in individual, resealable plastic carrying bags, complete with mask, reservoir bag, oxygen tube and pressure limiting valve 4.0 kPa H₂O.

The bags are color-coded for fast identification of infant, pediatric, and adult sizes.

Technical Specifications

Suitable for: baby 10-30 kg
 Stroke volume 1 hand 450 ml
 Resuscitator volume: 670 ml
 Size (length x diameter): 223x99 mm
 Bag reservoir volume: 2.600 ml