

# SILICONE RESUSCITATOR KIT with bag - adult



Product Code	34243
Unit of sale	1 pc
Minimum order	1
Type	Kit of medical devices
Class	

UK-REP	No
CH-REP	No

RDM (NSIS)	123113
EAN/UPC	8023279342437
GTIN	08023279342437

## Description

### SILICONE RESUSCITATOR BAG KIT - adult

Supplied with silicone face mask N° 4, 60 cm H<sub>2</sub>O pop-off valve, oxygen tube 2 metres, autoclavable reservoir valve and 2,500 ml reservoir bag (volume: 1,600 ml), helicoidal mouth wedge (34273), tongue holding forceps, 10 Guedel airways (34439) with Cubo emergency bag (27163).