

# SILICONE RESUSCITATOR BAG with MASK N 1 - infant



Product Code	34247
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2524894
CND	R03020201
EAN/UPC	8023279342475
UMDNS	13360
PARAF	903630440
GTIN	08023279342475

## Description

Transparent resuscitators made of silicone and polysulfone fully autoclavable at 134°C.

SILICONE RESUSCITATOR - infant

Supplied with silicone face mask N° 1 and 40 cm H<sub>2</sub>O pop-off valve. Volume: 280 ml.

Latex free.