


# PVC SINGLE USE RESUSCITATOR - infant with Pop-off valve



Product Code	34249
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2645407
CND	R03020201
EAN/UPC	8023279342499
UMDNS	13360
GTIN	08023279342499
GMDN	36086

## Description

PVC RESUSCITATOR - infant

Supplied with single use mask N° 1 and 40 cm H2O Pop-off valve - Volume: 280 ml

Disposable resuscitator bags with disposable mask, no return membrane valve and connection for PEEP valve.

PVC disposable, transparent mask shaped to allow a perfect adhesion to the patient's face. Latex free.