



| | |
|---------------|----------------|
| Product Code | 34255 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|----------------|
| RDM (NSIS) | 118947 |
| CND | R03010105 |
| EAN/UPC | 8023279342550 |
| UMDNS | 12447 |
| PARAF | 903630871 |
| GTIN | 08023279342550 |
| GMDN | 36066 |

Description

GIMA PLUS AUTOCLAVABLE FACE MASKS - ANAESTHESIA MASKS

Blue silicone inflating cuff and unbreakable transparent polysulfone hood.
They can be used with any resuscitator or lung ventilator.
Available in 5 sizes. Latex free.

AUTOCLAVABLE UP TO 120 TIMES.