



| | |
|---------------|----------------|
| Product Code | 34262 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|----|
| UK-REP | No |
|--------|----|

| | |
|---------|---------------|
| EAN/UPC | 8023279342628 |
| PARAF | 965618768 |
| GMDN | 17591 |

Description

SILICONE AUTOCLAVABLE RESUSCITATORS - WITH HANDLE

- Type: Infant
- Mask N. 1
- Pressure relief valve H2O: 40 cm
- Volume ml: 300

Transparent resuscitators made of medical grade silicone and fully autoclavable at 121°C.
Integrated handstrap for user comfort and complete control.

Made with no natural rubber or latex.

Multilingual manual: GB, FR, IT, ES, PT, SE, RO, GR.