

PU SHOWER STOOL with armrest



Product Code	42905
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

EAN/UPC	8023279429053
---------	---------------

Description

SHOWER STOOL WITH ARMRESTS

Shower stool made of silver aluminium frame with PU seat, foam armrest and adjustable foot tube.

Weight: 4.0 kg

Load: 136 kg

Size: 44.8-47.2x50-52.6x h 62.2-72 cm

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR.