

ROTATING SHOWER CHAIR with PU backrest and seat - load 136 kg



Product Code	42912
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

EAN/UPC	8023279429121
---------	---------------

Description

ROTATING SHOWER CHAIR

Height adjustable and 360° rotating seat chair.
Swing-up armrests make patient transfer easier. The lever locks the seat in 4 different positions.
Aluminium silver frame with comfortable PU padded seat and backrest.
Foam arm pad.
It can be folded.

Seat dimension: 40.7x33.7x h 5 cm
Max load: 136 kg.
Weight: 5.3 kg.
Packing size: 52.5x36x h 42.5 cm
Folding dimension: 51x34.5x h 40 cm

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, GR