

SHOWER COMMODE CHAIR



Product Code	42914
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC 6933257922860

Description

SHOWER COMMODE CHAIR - grey/white

Commode chair with double function: it can be used as a toilet and as a shower chair.
Made of aluminium and plastic material with PP frame for better load capacity.
Removable armrest and backrest, adjustable height in six levels.
It can be easily mounted without any tools.
Supplied with backrest and seat cushions, pan with lid and handle and foot pads.

Load capacity:
Weight: 5.1 kg.
Seat size: 46x43xh 42 cm.
Packing size: 55x47x24 cm.
Adjustable height