



Product Code	43045
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279430455
---------	---------------

PARAF	973923129
-------	-----------

Description

TRUSTY CANE - black

Sturdy, folding cane with built-in lights and extra wide pivoting head with 3 firm and skid-proof bases.

Height adjustable 85-97 cm.

Large base with Ø 9 cm.

Shaft in aluminium (thickness 1.2 mm)

Handle in ABS (15 cm).

Weight 420 g.

Need 3x1.5 V AAA batteries not included.

3 built-in LED lights, 1 forward and 2 down to light your way

Extra wide pivoting head revolves 360° degree

4 joints folding design

3 firm and skid-proof foot pads