

LARYNGOSCOPE SET 3 BL. MILLER 0-1, MC INT. 0 - paediatric



Product Code	34305
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2535677
CND	Z12021003
EAN/UPC	8023279343052
GTIN	08023279343052
GMDN	15076

Description

STANDARD RECHARGEABLE SETS (HANDLE + BLADES)

- Paediatric set (3 blades MILLER 0-1, Mc INTOSH 0)

Great features of the unique GIMA blades and handles

- standard 2.5 V vacuum bulb
- light ergonomic, satin finish handle to reduce glare
- Stainless Steel contacts guarantee long working life
- easy sterilization ETO or steam
- full range of Mc Intosh & Miller blades, latex free
- ISO 7376 fittings allow use of blades with existing handles
- set is packaged in an anti-shock case and in cardboard box
- 2.5 V handles work with both alkaline batteries or rechargeable batteries
- good light transmission of 3,500 LUX with 2.5 V handles
- autoclavable up to 134° for 5 minutes (approx. 2,000 times)
- non magnetic blades