



Product Code	43050
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279430509
PARAF	973923143

## Description

Innovative stick with safety light and alarm.

- front ultra-bright white 4 LED light up to 15 lumen, with adjustable positions
- ultra-bright red LED warning flash with 2 LED, up to 5 lumen
- ABS hand grip and aluminium frame
- 100 dB warning siren
- adjustable height: 10 positions from 78 to 100.5 cm

Works with 2 AAA batteries, not included.  
Maximum load: 110 kg

Supplied with carrying hand strap