



|               |                |
|---------------|----------------|
| Product Code  | 43072          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2184711       |
| CND        | Y120318       |
| EAN/UPC    | 8023279430721 |
| UMDNS      | 14393         |
| PARAF      | 977816406     |

## Description

QUADRIPOD Q3 with PU handle - bronze

Adjustable height: 70-88 cm - weight 0.5 kg

Maximum load: 100 kg