



Product Code	43076
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	15232
CND	Y120315
EAN/UPC	8023279430769
PARAF	987473055

Description

Walking aids for people who need to relieve the weight from the legs providing stability, higher than a traditional stick.

- TRIPOD WITH RING HANDGRIP - with adjustable height

Stick Structure: all the bases and the upper sliding tube are in light aluminium alloy. U-shaped tubes with strong injection-moulded plastic joint to provide strength and stiffness of the base.

Tube Locking clip: locking of the telescopic tubes with a clip equipped with a double pin; this improves the safety.

Tapered ring nut: supplied with a tapered ring nut allowing, by tightening, to eliminate noise.

Tips: high grip and have \varnothing of 35 mm, ISO EN 11334-4.

Ambidextrous ring shaped handgrip with non slip coating

Weight: 0.6 kg

Height: 79,5 to 93,5 cm

Max load: 100 kg

Made in Italy.