

# RUBBER TIPS tube 16 mm, base 35 mm



Product Code	43080
Unit of sale	box of 10 pcs
Minimum order	1
Type	Spare part
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279430806
PARAF	987470073

## Description

SPARE TIPS FOR STICKS, CRUTCHES, TRIPODS AND QUAD CANES

Tube: Ø 16 mm

Base: Ø 35 mm

Material: thermoplastic rubber

Wide range of spare tips that improve performance in term of comfort, grip, durability and adaptability to the ground.

Made in Italy.