

RUBBER TIPS tube 19 mm, base 55 mm



|               |              |
|---------------|--------------|
| Product Code  | 43086        |
| Unit of sale  | box of 2 pcs |
| Minimum order | 1            |
| Type          | Spare part   |
| Class         |              |

|        |    |
|--------|----|
| UK-REP | No |
| CH-REP | No |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279430868 |
| PARAF   | 987470123     |

Description

SPARE TIPS FOR STICKS, CRUTCHES, TRIPODS AND QUAD CANES

- Tube: Ø 19 mm  
Base: Ø 55 mm  
Material: pure rubber
- superior grip and high flexibility on irregular grounds
  - light shock absorbing effect

Wide range of spare tips that improve performance in term of comfort, grip, durability and adaptability to the ground.

Made in Italy.