



Product Code	43090
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	97984
CND	Y120306
EAN/UPC	3701127706949
PARAF	983513452

Description

Advance is the first polypropylene and elastomer bi-injected crutch, bringing enhanced comfort to the user.

Certified and tested up to 140 kg.

Adjustable in height 10 times per clip.

Integrated anti-noise cap and Pivoflex tip.

Made in France.

Technical Specifications

Weight kg: 0.51
 Height handle to floor cm: 77-98 cm
 Material of handle: polypropylene and elastomer
 Number of adjustment of height: 10
 Tube finishing/material: anodized aluminium
 Certified and tested for kg: 140
 Handgrip length cm: 10