



| | |
|---------------|----------------|
| Product Code | 43093 |
| Unit of sale | 1 pair |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|----|
| UK-REP | No |
|--------|----|

| | |
|---------|---------------|
| CND | Y120306 |
| EAN/UPC | 3701127751833 |
| PARAF | 951305186 |

Description

Advance is the first polypropylene and elastomer bi-injected crutch, bringing enhanced comfort to the user.

Certified and tested up to 140 kg.

Adjustable in height 10 times per clip.

Integrated anti-noise cap and Pivoflex tip.

Made in France.

Technical Specifications

Weight kg: 0.51
 Height handle to floor cm: 77-98 cm
 Material of handle: polypropylene and elastomer
 Number of adjustment of height: 10
 Tube finishing/material: anodized aluminium
 Certified and tested for kg: 140
 Handgrip length cm: 10