



Product Code	34311
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2837798
CND	Z12021099
EAN/UPC	8023279343113
PARAF	903147294
GTIN	08023279343113
GMDN	46827

Description

Blade N° 0 MILLER - premature/newborn
 Overall Length: 78 mm
 Distal Width: 11.5 mm

- Great features of the unique GIMA blades and handles
- standard 2.5 V vacuum bulb
 - light ergonomic, satin finish handle to reduce glare
 - Stainless Steel contacts guarantee long working life
 - easy sterilization ETO or steam
 - full range of Mc Intosh & Miller blades, latex free
 - ISO 7376 fittings allow use of blades with existing handles
 - set is packaged in an anti-shock case and in cardboard box
 - 2.5 V handles work with both alkaline batteries or rechargeable batteries
 - good light transmission of 3,500 LUX with 2.5 V handles
 - autoclavable up to 134° for 5 minutes (approx. 2,000 times)
 - non magnetic blades