



Product Code	43115
Unit of sale	1 pair
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes
RDM (NSIS)	97984
CND	Y120306
EAN/UPC	0370112770895
PARAF	975866207

Description

TIKI CRUTCHES - blue/orange - pair

For children or teenagers. The handle adjust from 16 to 21 cm.

Made in France.

Technical Specifications

Model TIKI
 Weight kg: 0.35
 Height handle to floor cm: 55-77
 Material of handle: Polypropylene
 Diameter upper/down pipe mm: 22/19
 Number of adjustment of height: 10 times adjustable every 20 mm
 Finishing/material: grey anodization/aluminium
 Certified and tested for kg: 100
 Standard used for testing: E N ISO 11334-1
 Adjustable handgrip: 3 times
 Handgrip size: Ø x l x h mm 32x100x160/210